



## Contents

<b>1</b>	<b>Introduction</b>	<b>2</b>
1.1	Basic Workflow Definition . . . . .	2
1.2	General Workflow Definition . . . . .	3
1.3	General Workflow - Complete Overview . . . . .	5
<b>2</b>	<b>Workflow - Single Steps</b>	<b>6</b>
2.1	Acquire Sample Data with Overview Scan . . . . .	6
2.2	Create Image Analysis Pipeline . . . . .	7
2.3	Define Overview Experiment . . . . .	9
2.4	Define Detailed Scan Experiment . . . . .	10
<b>3</b>	<b>Automation via OAD</b>	<b>11</b>
3.1	Pure Scripting or User Dialog . . . . .	11
3.2	User Dialog Requirements . . . . .	11
3.3	Prerequisites . . . . .	12
3.4	Create the User Dialog . . . . .	13
3.5	Overview Scan Experiment . . . . .	15
3.6	Find the interesting Objects . . . . .	16
3.7	Modify the Tile Dimensions . . . . .	18
3.8	Run the Detailed Scan . . . . .	19
<b>4</b>	<b>Testrun with FluoCells Slide</b>	<b>21</b>
<b>5</b>	<b>Test Run with Brain Slide</b>	<b>23</b>
<b>6</b>	<b>Appendix: Python Scripts</b>	<b>25</b>
6.1	Guided_Acquisition.czmac . . . . .	25
6.2	Custom DLLs . . . . .	30
<b>7</b>	<b>ToDo's and Limitations</b>	<b>30</b>
<b>8</b>	<b>Disclaimer</b>	<b>30</b>



## 1 Introduction

For a growing number of applications, it will be crucial to acquire data in a smart way. One way to achieve this goal is to build a smart microscope, which essentially means creating smart software workflows to control the hardware based on image analysis results.

Idea or Task:

- Scan or inspect a large area (or a long period of time).
- Detect an "interesting" object.
- Acquire detailed data for every event.
- Automate the workflow to minimize user interference.

First of all it is important to define what a **Object of Interest** can actually be.

- An object that meets specific criteria regarding its parameters (size, brightness, shape, intensity, ...)
- A specific change of a parameter during an experiment (e.g. a cell gets really "bright" upon ...)

It could be something quite simple. For instance one can have lots of cells, that are stained with blue dye, and only a few of them (maybe where the transfection worked ...) are also expressing GFP. The idea here would be to detect all cells that are positive for both colors and acquire an z-Stack for every cell (blue & green) that meets those criteria. Therefore this kind of application requires three major tasks:

1. Define the Overview Scan Experiment.
2. Define the object detection rules, e.g. setup image analysis.
3. Define the Detailed Scan(s) to be carried out in case of a "positive" object.

### 1.1 Basic Workflow Definition

The main workflow can be summarized as follows:

1. Acquire some sample data showing a object of interest.
2. Setup an Image Analysis Pipeline to detect those events.
3. Define an experiment which does the Overview Scan.
4. Define an experiment which does the Detailed Scan.
5. Automate the entire workflow.

The goal of this tutorial is to create an automated workflow that can be used to easily setup a **Guided Acquisition**. This requires some knowledge about the OAD macro environment and its scripting language Python.

## 1.2 General Workflow Definition

This is the outline of the general workflow, which would be automated.




Figure 1: General workflow for Guided Acquisition.

A different way to visualize this is shown in the figure below.




Figure 2: Larger area with some candidates for a "interesting" objects.

### 1.3 General Workflow - Complete Overview




Figure 3: Workflow - Actionable Information 1.




Figure 4: Workflow - Actionable Information 2



## 2 Workflow - Single Steps

This part of the tutorial will explain all required steps to set up the complete workflow in more detail. Some knowledge about image analysis is required here.

### 2.1 Acquire Sample Data with Overview Scan

The first crucial step is to have a basic idea of what the actual rare event will look like, inside an real image acquired, with the appropriate parameters (light intensity, detector settings, filters, objective, ...). An ideal sample data set should be acquired using the "real" acquisition parameters will be used to define the overview scan experiment later on.

Typically such an overview scan is acquired using a lower magnification in combination with a tile experiment. The exact parameters will be specific for the application, but the main idea is always the same:

**The overview image must contain the information to locate the objects of interest based on "some" features that can be retrieved via an appropriate image analysis.**

Once a representative sample data set is available, one can start setting up an image analysis pipeline. ZEN offers currently two ways do to so:

1. The easy way – Use the Analysis Wizard.
2. Program your own image analysis pipeline using an OAD macro (this is not covered by this tutorial).

The most comfortable way is to use the wizard; this offers enough flexibility to cover most of the applications.

## 2.2 Create Image Analysis Pipeline

Here one already sees a "half-ready" image analysis pipeline. In order to automate, the wizard step containing the feature selection is crucial.

In order to relocate all detected objects, it is of course required to retrieve the current XY(Z) position of every detected event.

The main idea is to determine all parameters required to get access to the stage coordinates of a detected object.




Figure 5: Select the required positional features inside the wizard.



Currently one needs to select the following parameters:

- Object ID
- Bound Center Stage X and Y
- Bound Width and Height

The parameters **Bound Center Stage X** and **Bound Stage Center Y** yield absolute stage coordinates of the detected objects. These will be used to relocate the objects for the detailed scan later on.

The parameters **Bound Width** and **Bound Height** yield in the absolute width and height of the bounding rectangle. This is required if the detected object is bigger than a single frame. In such a case the tile region of the detailed experiment will be adapted to the size of the bounding rectangle automatically.

Depending on your application, it might be useful to measure additional parameters. Feel free to add whatever might be useful.

## 2.3 Define Overview Experiment

Usually this is done using the Tiles and Positions module to scan a large area. This tutorial assumes that one is already familiar with this ZEN module.

**Important Remark:** Since one wants to relocate the detected objects, it is required to calibrate the XY stage before the overview scan, and a rigid and stiff sample holder should be used.

For this tutorial one can use an already acquired image representing the real overview scan (instead of a real experiment). The shown image is a 16x13 tile image acquired with a 10X magnification. The result of such a scan followed by the image analysis is shown here:




Figure 6: Result of the image analysis on the overview scan image.

## 2.4 Define Detailed Scan Experiment

First of all it is important to understand, that there is no such thing as a typical "Detailed Scan". What this experiment (or even a workflow) might be, depends on the application. In a general sense such a detailed scan is "some kind" of experiment carried at a specific position based upon the results of the image analysis. Examples could be:

- Simple Z-Stack with a high-NA objective lens.
- Multi-Channel Z-Stack using an optical sectioning method like SD, Apotome or AiryScan
- One of the above combined with a tile experiment.
- A series of subsequent experiments with different image modalities.

Since this concept is based on OAD, it is possible to automate almost any kind of workflow at a given position.




Figure 7: Setup of a simple detailed scan experiment.

Now the preparation is complete. So far we have available:

- The image analysis pipeline based on some sample data
- The Overview Scan Experiment = tile experiment created by the Tiles & Positions Module
- The Detailed Scan Experiment, which is also a tile experiment.

And now the real interesting part begins – how to setup a complete workflow out of those building blocks using an OAD Python script inside ZEN Blue.



## 3 Automation via OAD

### 3.1 Pure Scripting or User Dialog

First of all one has to decide if a pure script is sufficient or if there should be some kind of simple user interface (very basic wizard).

- A pure script is of course the faster solution, but is may be difficult to use for less advanced users.
- A basic wizard offers less flexibility, but may be easier to use for most users.

The choice clearly depends on the user, since there is no right or wrong approach here. For this tutorial we will focus on the 2nd approach. Once one has figured out how to realize this, the 1st approach will be straight forward since the basic workflow is identical.

### 3.2 User Dialog Requirements

The basic user interface of the dialog should have the following functionality:

- **Select the overview scan experiment.**
- **Select the image analysis pipeline** (to analyze the overview image)
- **Select the detailed scan experiment** (to acquire at "found" positions)
- **Define options for focussing** (depends on available hardware)
- **Select a folder to store the image data and analysis results.**



### 3.3 Prerequisites

Now we are ready to create the OAD macro required to control the complete workflow. The first step is (as always) to import the required modules and define the folder locations, etc. Additionally a python function is defined here for later usage. It will check, if a chosen directory is really empty before the actual workflow starts.

```
21 # import custom DLLs required for Rare Event Detection
22 import clr
23 clr.AddReference('RETools.dll')
24 clr.AddReference('TileTools.dll')
25 import RETools
26 import Tiles
27
28 from System.IO import File, Directory, Path
29 import sys
30
31 # optional debugging output
32 verbose = True
33 # version number for dialog window
34 version = 4.1
35 # file name for overview scan
36 ovscan_name = 'OverviewScan.czi'
37 # additional XY offset for possible LSM port relative to the Camera port (zero)
38 dx_LSM = 0.0
39 dy_LSM = 0.0
40
41 def dircheck(folder2check, createdir=False):
42
43     # check if the destination folder is empty
44     direxists = Directory.Exists(folder2check)
45     if direxists:
46         isEmpty = Directory.GetFiles(folder2check).Length
47         # check if the directory is empty in case it already exists
48         if isEmpty != 0:
49             print 'Directory already contains files! Confirm Selection or Create new Directory.'
50             SelectFolderDialog = ZenWindow()
51             SelectFolderDialog.Initialize('New Folder Selection? Directory Contains Files!', 650, 100, True, True)
52             # add components to dialog
53             SelectFolderDialog.AddFolderBrowser('fd', 'Confirm Folder or Select new Directory', folder2check, '0', '0')
54             sf = SelectFolderDialog.Show()
55             if sf.HasValue:
56                 message = 'Macro was canceled by user.'
57                 print message
58                 raise SystemExit
59
60             folder2check = str(sf.GetValue('fd'))
61             print 'Creating Output Folder: ', folder2check
62
63     # create directory when option is set
64     if createdir:
65         print 'Creating Output Folder: ', folder2check
66         Directory.CreateDirectory(folder2check)
67
68     return folder2check
69
70 # clear console output
71 Zen.Application.MacroEditor.ClearMessages()
72
73 # check the location of folder where experiment setups and image analysis settings are stored
74 docfolder = Zen.Application.Environment.GetFolderPath(ZenSpecialFolder.UserDocuments)
75 #imgfolder = Zen.Application.Environment.GetFolderPath(ZenSpecialFolder.ImageAutoSave)
76 imgfolder = r'D:\RareEvent_FL-Slide1\Output_RareEvent'
77
78 # get list with all existing experiments and image analysis setup
79 expfiles = RETools.RareEventTools.GetExperimentSetups(docfolder)
80 ipfiles = RETools.RareEventTools.GetImageAnalysisSetups(docfolder)
```



### 3.4 Create the User Dialog

The next task is to create the user dialog using the requirements referred to earlier on. Once the dialog is closed, the results have to be collected.

Additionally one must check if the Detailed Scan experiment is a tile experiment, where the tile size has to be adapted accordingly. This is done using a little tool function.

```
83 # collect information about the configured objectives and optovars
84 objectives = RETools.RareEventTools.CreateObjectiveInfo()
85 optovars = RETools.RareEventTools.CreateOptovarInfo()
86
87 # Activate Zen.Application
88 RareEventDialog = ZenWindow()
89 RareEventDialog.Initialize('Guided Acquisition - Version : ' + str(version), 650, 750, True, True)
90 # add components to dialog
91 RareEventDialog.AddLabel('1) Select Overview Experiment -----')
92 RareEventDialog.AddDropDown('overview_exp', 'Overview Scan Experiment', expfiles, 0)
93 RareEventDialog.AddDropDown('objectiveOV', 'Objective Overview Scan', objectives.Names, 1)
94 RareEventDialog.AddDropDown('optovarOV', 'After-Mag OverViewScan', optovars.Names, 2)
95 RareEventDialog.AddCheckbox('SWAF_before_overview', 'OPTION - (Find Surface) & SWAF before Overview', True)
96 RareEventDialog.AddLabel('2) Select Image Analysis to detect objects -----')
97 RareEventDialog.AddDropDown('ip_pipe', 'Image Analysis Pipeline', ipfiles, 0)
98 RareEventDialog.AddLabel('3) Select Detail Scan Experiment -----')
99 RareEventDialog.AddDropDown('detailed_exp', 'Detailed Scan Experiment', expfiles, 1)
100 RareEventDialog.AddDropDown('objectiveDT', 'Objective Detailed Scan', objectives.Names, 2)
101 RareEventDialog.AddDropDown('optovardt', 'After-Mag Detailed Scan', optovars.Names, 2)
102 RareEventDialog.AddCheckbox('checkoffset', 'OPTION - Check offset between FS and SWAF', False)
103 RareEventDialog.AddCheckbox('manualoffset', 'OPTION - Enter Offset manually', False)
104 RareEventDialog.AddDoubleRange('offsetvalue', 'Enter manual Offset [micron]', 0, 0, 100)
105 RareEventDialog.AddLabel('4) Specify Output Folder to save the images -----')
106 RareEventDialog.AddFolderBrowser('outfolder', 'Output Folder for Images and Data Tables', imgfolder)
107
108 # show the window
109 result = RareEventDialog.Show()
110 if result.HasCanceled:
111     message = 'Macro was canceled by user.'
112     print message
113     raise SystemExit
114
115 # get the values and store them
116 OverViewExpName = str(result.GetValue('overview_exp'))
117 ImageAS = str(result.GetValue('ip_pipe'))
118 DetailExpName = str(result.GetValue('detailed_exp'))
119 OutputFolder = str(result.GetValue('outfolder'))
120 SWAF_beforeOV = result.GetValue('SWAF_before_overview')
121 CheckOffset = result.GetValue('checkoffset')
122 ManualOffset = result.GetValue('manualoffset')
123 OffsetValue = result.GetValue('offsetvalue')
124 objOV = result.GetValue('objectiveOV')
125 optoOV = result.GetValue('optovarOV')
126 objDT = result.GetValue('objectiveDT')
127 optoDT = result.GetValue('optovardt')
128
129 # print values - this is optional
130 print 'Overview Scan Experiment : ' + OverViewExpName
131 print 'Combination OverView Scan: ', objOV + ' + ' + optoOV
132 print 'Image Analysis Pipeline : ' + ImageAS
133 print 'Detailed Scan Experiment : ' + DetailExpName
134 print 'Combination Detailed Scan: ', objDT + ' + ' + optoDT
135 print 'Output Folder for Data : ' + OutputFolder, '\n'
136
137 # check directory
138 OutputFolder = dircheck(OutputFolder, createdir=True)
```

The resulting dialog window is shown below.



Figure 8: Simple User Dialog to run the Guided Acquisition WorkFlow

The dialog has four main steps where the user must define the desired options:

1. Select **overview experiment** and additional focus options.
2. Select **image analysis** to detect the objects.
3. Select **detailed experiment** and additional focus options.
4. Choose an **output folder** to save the data.



### 3.5 Overview Scan Experiment

At this point we have all the information we need to actually start the workflow. The first step is to execute the actual overview scan experiment, which is most likely a tile experiment. The resulting image will be saved to the specified folder.

```
157 ##### START OVERVIEW SCAN EXPERIMENT #####
158
159 # use the correct objective and optovar for the overview scan
160 RETools.RareEventTools.SetObjectiveByName(objOV, False)
161 RETools.RareEventTools.SetOptovarByName(optoOV, False)
162
163 if SWAF_beforeOV==True:
164     try:
165         # initial focussing via FindSurface to assure a good starting position
166         # requires DF2 for Observer or Celldiscoverer ? --&gt; otherwise comment line !!!
167         Zen.Acquisition.FindSurface()
168     except:
169         print 'Was not able to run Find Surface.'
170
171     try:
172         # run the SWAF using the settings from the OVScan --&gt; check SWAF for overview experiment !!!
173         Zen.Acquisition.FindAutofocus(OVScan)
174     except:
175         print 'Was not able to run SWAF using the seetings: ', OverViewExpName
176
177 # get the resulting z-position
178 znew = Zen.Devices.Focus.ActualPosition
179
180 # try to adapt the Tile Experiment with new Z-Position
181 try:
182     Tiles.TileTools.ModifyTileRegionsZonly(OVScan, znew)
183     print 'Adapted Z-Position of Tile OverView. New Z = ',znew
184 except:
185     print 'Was not able to adapt Z-Position of Overview Scan Experiment.'
186
187 # execute the experiment
188 print 'Running OverviewScan Experiment.\n'
189 output_OVScan = Zen.Acquisition.Execute(OVScan)
190 OVScan.Close()
191 # For testing purposes - Load overview scan image automatically instead of executing the "real" experiment
192 #output_OVScan =
193 #    Zen.Application.LoadImage(r'c:\Users\MI1SRH\Documents\Testdata_Zeiss\RareEvent_Test_Wizard\OverViewScan_Test_raw.czi',
194 #    False)
195
196 # show the overview scan inside the document area
197 Zen.Application.Documents.Add(output_OVScan)
198 ovdoc = Zen.Application.Documents.GetByName(output_OVScan.Name)
199
200 # save the overview scan image inside the select folder
201 output_OVScan.Save(OutputFolder + '\\\\' + ovscan_name)
202 ###### END OVERVIEW SCAN EXPERIMENT #####
```

A very important part of such an automated workflow is to find the focus. To ensure this, the desired objective and optovar are changed before the overview starts, followed by a series of **optional** focus actions:

1. Find the surface of the sample carrier.
2. Focus onto the specimen sample via software autofocus.
3. Store the resulting new Z-position.
4. Modify the tile overview experiment using the new Z-position.

When testing out applications like this it is very useful to use test data sets to save time instead of re-running an acquisition many times. Therefore it can be helpful to **comment** the lines where the real experiment is executed and **uncomment** the line inside the script, where one can load an already acquired overview image.



### 3.6 Find the interesting Objects

Now it is time to run the image analysis on the overview image. As a result, two tables will be created inside ZEN. The first contains information about all objects and the second contains information about the single objects (for instance their stage positions). This second table is what will be used as an input for the detailed scan later on.

```
203 # Load analysis setting created by the wizard or an separate macro
204 ias = ZenImageAnalysisSetting()
205 # for simulation use: 000 - RareEventExample.czias
206 ias.Load(ImageAS)
207 # Analyse the image
208 Zen.Analyzing.Analyze(output_OVScan,ias)
209 # Create Zen table with results for all detected objects (parent class)
210 A11Obj = Zen.Analyzing.CreateRegionsTable(output_OVScan)
211 # Create Zen table with results for each single object
212 SingleObj = Zen.Analyzing.CreateRegionTable(output_OVScan)
213
214 # check for existence of required column names
215 DetailScan = Zen.Acquisition.Experiments.GetByName(DetailExpName)
216 DetailScan.SetActive()
217 DetailScanExp = Tiles.TileTools.IsTilesExperiment(DetailScan)
218 out = RETools.RareEventTools.CheckColumns(SingleObj, DetailIsTileExp)
219 DetailScan.Close()
220
221 # 1st item is a bool indicating if all required columns could be found
222 columnsOK = out.Item1
223
224 if columnsOK == False:
225     print 'Execution stopped. Required Columns are missing.'
226     raise Exception('Execution stopped. Required Columns are missing.')
227
228 # 2nd item is a dictionary containing the tile properties
229 # colID: used keys = 'BCcolx' / 'BCcoly' / 'BCwidthcolx' / 'BCheightcoly'
230 colID = out.Item2
231
232 # show and save data tables to the specified folder
233 Zen.Application.Documents.Add(A11Obj)
234 Zen.Application.Documents.Add(SingleObj)
235 A11Obj.Save(OutputFolder + '\\\\OverviewTable.csv')
236 SingleObj.Save(OutputFolder + '\\\\SingleObjectsTable.csv')
237
238 # check the number of detected objects = rows inside image analysis table
239 num_POI = SingleObj.RowCount
240
241 # switch to new magnification and optovar before the detailed scan starts
242 RETools.RareEventTools.SetObjectivebyName(objDT, False)
243 RETools.RareEventTools.SetOptovarbyName(optoDT, False)
244
245 # determine offset between FindSurface and objects found by SWAF -->; required DF2
246 if CheckOffset==True:
247     # move to the first valid XY position to determine the offset
248     xpos_1st = SingleObj.GetValue(0, colID['BCcolx']) # get x stage position from list
249     ypos_1st = SingleObj.GetValue(0, colID['BCcoly']) # get y stage position from list
250     Zen.Devices.Stage.MoveTo(xpos_1st + dx_LSM, ypos_1st + dy_LSM)
251
252     # get the actual offset using FindSurface followed by SWAF -->; requires DF2
253     try:
254         dzFS = RETools.RareEventTools.getOffsetFindSurfaceSWAF(DetailExpName, False)
255         print 'Automatic Offset dzFS measured [micron]: ', round(dzFS, 2)
256         # store this new value inside the DF in order to use RecallFocus later on
257         print 'Store Z-Value Offset for RecallFocus.\n'
258         # store offset inside DF to use it with Recall Focus using DF2
259         Zen.Acquisition.StoreFocus()
260         useRecallFocus = True
261     except:
262         print 'Automatic Offset Check failed. Setting Offset dzFS = 0.'
263         dzFS = 0.0
264         useRecallFocus = False
265
266     elif CheckOffset==False:
267         dzFS = 0.0
268         useRecallFocus = False
269
270 if ManualOffset==True and CheckOffset==False:
271     # use the manually entered offset only when dz was not measured automatically
272     dzFS = OffsetValue
273     print 'Manual Offset [micron]: ', round(dzFS, 2)
```

As shown in figure 5, it is required to have the correct image analysis features selected. Only the will the respective columns inside the table exist. This has to be checked because otherwise the workflow will be impossible.

The "checking" itself is again done by a tool function. The output will be a boolean and a dictionary in order to look up the column index for the respective image analysis parameters. This has the advantage of being independent from a specific column order. One just has to make sure that all required columns are there.

If the option inside the initial dialog window was checked, the offset between the sample carrier surface and the actual specimen will be measured and stored inside the Definite Focus via **Store Focus**, so that **Recall Focus** can be used later on when needed. The dialog also offers the possibility to enter an arbitrary offset manually, which will only be applied if the offset was not determined automatically.




Figure 9: The results of the overview scan and the image analysis.

Finally the number of detected objects is checked by looking up the number of rows inside the table containing the information about the single objects.



### 3.7 Modify the Tile Dimensions

Once one has the list of objects with the stage coordinates, the next task is to create a loop to cycle through all detected positions.

- Move to the correct XYZ stage positions.
- Run Focus Actions when required.
- Initialize the DetailScan experiment at the new position.
- Adapt TilePosition using *Tiles.TileTools.ModifyTileRegionsXYZ*.
- Modify TileSize using *Tiles.TileTools.ModifyTileRegionsSize*.

```
278 ##### START DETAILED SCAN EXPERIMENT #####
279
280 # execute detailed experiment at the position of every detected object
281 for i in range(0, num_POI, 1):
282
283     # get the object information from the position table
284     POI_ID = SingleObj.GetValue(i, 0) # get the ID of the object - IDs start with 2 !!!
285     xpos = SingleObj.GetValue(i, colID['BCcolx']) # get X-stage position from table
286     ypos = SingleObj.GetValue(i, colID['BCcoly']) # get Y-stage position from table
287
288     # move to the current position
289     Zen.Devices.Stage.MoveTo(xpos + dx_LSM, ypos + dy_LSM) # comment this, if one uses a simulated experiment
290     print 'Moving Stage to Object ID:', POI_ID, ' at :', round(xpos, 2), round(ypos, 2)
291
292     if useRecallFocus == False:
293         # Initial FindSurface before the Detail Scan starts
294         try:
295             Zen.Acquisition.FindSurface()
296         except:
297             'Was not able to run Find Surface.'
298             # calculate new focus position plus offset and move z-drive
299             zpos = Zen.Devices.Focus.ActualPosition + dzFS
300             print 'Move to Z-Position: ', round(zpos, 2)
301             Zen.Devices.Focus.MoveTo(zpos)
302
303     elif useRecallFocus == True:
304         # alternative solution - use RecallFocus
305         try:
306             Zen.Acquisition.RecallFocus()
307         except:
308             print 'Was not able to run Recall Focus.'
309             zpos = Zen.Devices.Focus.ActualPosition
310             print 'New z-position after Recall Focus: ', zpos
311
312     # load the predefined detailed scan experiment
313     DetailScan = Zen.Acquisition.Experiments.GetByName(DetailExpName)
314
315     # only modify the Tile Properties if required IAS features BoundWidth and BoundHeight were found
316     # if experiment is a Tile Experiment
317     if DetailIsTileExp == True:
318
319         # Modify tile center position - get bounding rectangle width & height in microns
320         bcwidth = SingleObj.GetValue(i, colID['BCWidthcolx'])
321         bcheight = SingleObj.GetValue(i, colID['BCHeightcoly'])
322         print 'Width and Height : ', str(round(bcwidth, 2)), str(round(bcheight, 2))
323         print 'Modifying Tile Properties XYZ Position and width & height.'
324         # Modify the XYZ position of the tile region on-the-fly
325         Tiles.TileTools.ModifyTileRegionsXYZ(DetailScan, xpos, ypos, zpos)
326         # Modify ConturSize for the tile according to the size of the bounding rectangle
327         Tiles.TileTools.ModifyTileRegionsSize(DetailScan, bcwidth, bcheight)
328         print 'New Tile Properties: ', round(xpos, 2), round(ypos, 2), round(zpos, 2), round(bcwidth, 2), round(bcheight, 2)
```

### 3.8 Run the Detailed Scan

The last steps, which need to be executed now, are:

- Run the (modified) detailed scan experiment here.
- Save the data to the specified folder.
- Close the image document in ZEN.
- Rename the file using the current object ID.

```

330 # execute the experiment
331 print 'Running Detail Scan Experiment at new XYZ position.'
332

```




Figure 10: The results of the detailed scan - tile regions adapts to objects.

If the optional output(verbose = True) option was used, one can see what is going on inside the IDE Watch viewblock.

```

Overview Scan Experiment : 000 - RareEvent_OverviewScan.czexp
Image Analysis Pipeline : 000 - RareEventExample.czias
Detailed Scan Experiment : 001 - RareEvent_DetailedScan.czexp
Output Folder for Data : C:\MyData\RareEvent_Output

Detailed Experiment is a Tile Experiment.

Running OverviewScan Experiment.

Detected Column ID at Col: 1
Detected Column BoundCenterXStage at Col: 2
Detected Column BoundCenterYStage at Col: 3
Detected Column BoundWidth at Col: 4
Detected Column BoundHeight at Col: 5
Detected Object ID : 2 XY Position : 85450.18 -17317.87
Width & Height : 918.48 722.4
Modifying Tile Properties...
Tile Center Position XYZ : 85450.178 -17317.866 0.0
Tile Contour Size : 1815.472 1089.366
Tile Size XY : 8 x 5

```

```

Tile Center Position XYZ : 85450.178 -17317.866 0.0
Tile Contour Size : 918.48 722.4
Tile Size XY : 4 x 4
Moving Stage to Object ID: 2 at : 85450.178 -17317.866
Running DetailedExperiment at new position.

Renaming File: Experiment-52.czi to: DTScan_ID2.czi
Detected Object ID : 3 XY Position : 84258.22 -16610.95
Width & Height : 784.32 732.72
Modifying Tile Properties...
Tile Center Position XYZ : 84258.218 -16610.946 0.0
Tile Contour Size : 918.48 722.4
Tile Size XY : 4 x 4
Tile Center Position XYZ : 84258.218 -16610.946 0.0
Tile Contour Size : 784.32 732.72
Tile Size XY : 4 x 4
Moving Stage to Object ID: 3 at : 84258.218 -16610.946
Running DetailedExperiment at new position.

Renaming File: Experiment-53.czi to: DTScan_ID3.czi
Detected Object ID : 4 XY Position : 82756.66 -15578.95
Width & Height : 856.56 505.68
Modifying Tile Properties...
Tile Center Position XYZ : 82756.658 -15578.946 0.0
Tile Contour Size : 784.32 732.72
Tile Size XY : 4 x 4
Tile Center Position XYZ : 82756.658 -15578.946 0.0
Tile Contour Size : 856.56 505.68
Tile Size XY : 4 x 3
Moving Stage to Object ID: 4 at : 82756.658 -15578.946
Running DetailedExperiment at new position.

Renaming File: Experiment-54.czi to: DTScan_ID4.czi
All Positions done. RareEventDetection finished.

```

That is it. All image data sets acquired at the detected positions are stored inside the correct folder using the object IDs as names. The tables are saved in the same location.




Figure 11: The results of the detailed scan stored inside a folder.

## 4 Testrun with FluoCells Slide

To test the workflow on a real sample one can use the FluoCells slide. The general idea is to setup an image analysis that detects some "interesting" objects. In order to create a somewhat short test run it is recommended to use an image analysis pipeline that only yields in a few positive detection events.




Figure 12: The results of the overview scan and image analysis.




Figure 13: Detailed view on a detected object.




Figure 14: Detailed images from some of the detected objects.

The image analysis yielded **24** positive detections of cells with a large nucleus, preferably cells which were "caught" during cell division. Keep in mind, that the image analysis was not "fine-tuned" to catch all cell divisions. It is all about testing out the general workflow. For this reason it is also good practice to acquire a small overview scan at the beginning, e.g. a 2x2 tile image.

## 5 Test Run with Brain Slide

The same basic workflow can also be used to detect brain slices on a slide. All one needs to change is the image analysis pipeline in order to detect the slices.




Figure 15: The results of the overview scan and image analysis.




Figure 16: A more detailed image of the first specimen with three zoomed-in regions.



The image analysis yielded **eight** positive regions for brains slices on this slide. The properties of the bounding rectangle were used to modify the required tile experiment. The final detailed scan contained 36 tiles and it was imaged using only a hardware-based focus system.

Keep in mind that the actual focusing during such a experiment must be still part of the focus strategies with the Detail Scan experiment. Which one must be used here greatly depends on the nature of the sample and the actual application. The Guided Acquisition tool only supports finding the initial positions based on the overview scan image.



## 6 Appendix: Python Scripts

### 6.1 Guided\_Acquisition.czmac

This is the complete ZEN Blue python script used to realize the workflow described above. In principle it can work on any motorized microscope stand running under ZEN blue (with some modifications), but it was tested and optimized mainly for the Celldiscoverer 7 and for the AxioObserver 7. Both of those stands have the Definite Focus 2 as an additional hardware option, which is strongly advised for such highly automated workflows, where fast and efficient focus options, such as Find Surface are really a big plus.

```
1  """
2  Author: Sebastian Rhode
3  Date: 2016_11_29
4  File: Guided_Acquisition_DLLs.czmac
5  Version: 4.2
6
7  Optimized for the use with Celldiscoverer 7 and DF2.
8  Please adapt Focussing commands, especially FindSurface when using with other stands
9
10 1) - Select Overview Scan Experiment
11 2) - Select appropriate Image Analysis Pipeline
12 3) - Select Detailed Scan Experiment
13 4) - Specify the output folder for the image and data tables
14
15 Requires the following DLLs to be found inside the Zen program folder:
16
17 - RETools.dll
18 - TileTools.dll
19 """
20
21 # import DLLs
22 import clr
23 clr.AddReference('RETools.dll')
24 clr.AddReference('TileTools.dll')
25 import RETools
26 import Tiles
27 import time
28
29 from System.IO import File, Directory, Path
30 import sys
31
32 # optional debugging output
33 verbose = True
34 # version number for dialog window
35 version = 4.2
36 # file name for overview scan
37 ovscan_name = 'OverviewScan.czii'
38 # additional XY offset for possible LSM port relative to the Camera port (zero)
39 dx_LSM = 0.0
40 dy_LSM = 0.0
41
42 def dircheck(folder2check, createdir=False):
43
44     # check if the destination folder is empty
45     direxists = Directory.Exists(folder2check)
46     if direxists:
47         isEmpty = Directory.GetFiles(folder2check).Length
48         # check if the directory is empty in case it already exists
49         if isEmpty != 0:
50             print 'Directory already contains files! Confirm Selection or Create new Directory.'
51             SelectFolderDialog = ZenWindow()
52             SelectFolderDialog.Initialize('New Folder Selection? Directory Contains Files!', 650, 100, True, True)
53             # add components to dialog
54             SelectFolderDialog.AddFolderBrowser('fd', 'Confirm Folder or Select new Directory', folder2check)
55             sf = SelectFolderDialog.Show()
56             if sf.HasCanceled:
57                 message = 'Macro was canceled by user.'
58                 print message
59                 raise SystemExit
60
61             folder2check = str(sf.GetValue('fd'))
62             print 'Creating Output Folder: ', folder2check
```



```
63      # create directory when option is set
64      if createdir:
65          print 'Creating Output Folder: ', folder2check
66          Directory.CreateDirectory(folder2check)
67
68      return folder2check
69
70
71      # clear console output
72      Zen.Application.MacroEditor.ClearMessages()
73
74      # check the location of folder where experiment setups and image analysis settings are stored
75      docfolder = Zen.Application.Environment.GetFolderPath(ZenSpecialFolder.UserDocuments)
76      #imgfolder = Zen.Application.Environment.GetFolderPath(ZenSpecialFolder.ImageAutoSave)
77      imgfolder = r'c:\Output\Guided_Acquisition'
78
79      # get list with all existing experiments and image analysis setup
80      expfiles = RETools.RareEventTools.GetExperimentSetups(docfolder)
81      ipfiles = RETools.RareEventTools.GetImageAnalysisSetups(docfolder)
82
83      # collect information about the configured objectives and optovars
84      objectives = RETools.RareEventTools.CreateObjectiveInfo()
85      optovars = RETools.RareEventTools.CreateOptovarInfo()
86
87      # Activate Zen.Application
88      RareEventDialog = Zenwindow()
89      RareEventDialog.Initialize('Guided Acquisition - Version : ' + str(version), 650, 750, True, True)
90      # add components to dialog
91      RareEventDialog.AddLabel('1) Select Overview Experiment -----')
92      RareEventDialog.AddDropDown('overview_exp', 'Overview Scan Experiment', expfiles, 0)
93      RareEventDialog.AddDropDown('objectiveOV', 'Objective Overview Scan', objectives.Names, 1)
94      RareEventDialog.AddDropDown('optovarOV', 'After-Mag OverViewScan', optovars.Names, 2)
95      RareEventDialog.AddCheckbox('SWAF_before_overview', 'OPTION - (Find Surface) & SWAF before Overview', True)
96      RareEventDialog.AddLabel('2) Select Image Analysis to detect objects -----')
97      RareEventDialog.AddDropDown('ip_pipe', 'Image Analysis Pipeline', ipfiles, 0)
98      RareEventDialog.AddLabel('3) Select Detail Scan Experiment -----')
99      RareEventDialog.AddDropDown('detailed_exp', 'Detailed Scan Experiment', expfiles, 1)
100     RareEventDialog.AddDropDown('objectiveDT', 'Objective Detailed Scan', objectives.Names, 2)
101     RareEventDialog.AddDropDown('optovarDT', 'After-Mag Detailed Scan', optovars.Names, 2)
102     RareEventDialog.AddCheckbox('checkboxoffset', 'OPTION - Check offset between FS and SWAF', False)
103     RareEventDialog.AddCheckbox('manualoffset', 'OPTION - Enter Offset manually', False)
104     RareEventDialog.AddDoubleRange('offsetvalue', 'Enter manual Offset [micron]', 0, 0, 100)
105     RareEventDialog.AddLabel('4) Specify Output Folder to save the images -----')
106     RareEventDialog.AddFolderBrowser('outfolder', 'Output Folder for Images and Data Tables', imgfolder)
107
108     # show the window
109     result = RareEventDialog.Show()
110     if result.HasCanceled:
111         message = 'Macro was canceled by user.'
112         print message
113         raise SystemExit
114
115     # get the values and store them
116     OverViewExpName = str(result.GetValue('overview_exp'))
117     ImageAS = str(result.GetValue('ip_pipe'))
118     DetailExpName = str(result.GetValue('detailed_exp'))
119     OutputFolder = str(result.GetValue('outfolder'))
120     SWAF_beforeOV = result.GetValue('SWAF_before_overview')
121     CheckOffset = result.GetValue('checkboxoffset')
122     ManualOffset = result.GetValue('manualoffset')
123     OffSetValue = result.GetValue('offsetvalue')
124     objOV = result.GetValue('objectiveOV')
125     optoOV = result.GetValue('optovarOV')
126     objDT = result.GetValue('objectiveDT')
127     optoDT = result.GetValue('optovarDT')
128
129     # print values - this is optional
130     print 'Overview Scan Experiment : ' + OverViewExpName
131     print 'Combination OverView Scan: ', objOV + ' + ' + optoOV
132     print 'Image Analysis Pipeline : ' + ImageAS
133     print 'Detailed Scan Experiment : ' + DetailExpName
134     print 'Combination Detailed Scan: ', objDT + ' + ' + optoDT
135     print 'Output Folder for Data : ' + OutputFolder, '\n'
136
137     # check directory
138     OutputFolder = dircheck(OutputFolder, createdir=True)
139
140     # check Detail Scan experiment for tiles
141     DetailScan = Zen.Acquisition.Experiments.GetByName(DetailExpName)
142     DetailScan.SetActive()
143     DetailIsTileExp = Tiles.TileTools.IsTilesExperiment(DetailScan)
144     DetailScan.Close()
145
146     # check Overview Scan for tiles
147     OVScan = Zen.Acquisition.Experiments.GetByName(OverViewExpName)
148     OVScan.SetActive()
149     OVScanIsTileExp = Tiles.TileTools.IsTilesExperiment(OVScan)
150
151     if (OVScanIsTileExp == False or DetailIsTileExp == False):
```



```
152     message = 'Overview or Detail Scan are not a Tile Experiment.'
153     print message
154     raise SystemExit
155     OVScan.Close()
156
157 ##### START OVERVIEW SCAN EXPERIMENT #####
158
159 # use the correct objective and optovar for the overview scan
160 RETools.RareEventTools.SetObjectiveByName(objOV, False)
161 RETools.RareEventTools.SetOptovarByName(optoOV, False)
162
163 if SWAF_beforeOV==True:
164     try:
165         # initial focussing via FindSurface to assure a good starting position
166         # requires DF2 for Observer or Celldiscoverer 7 --&gt; otherwise comment line !!!
167         Zen.Aquisition.FindSurface()
168     except:
169         print 'Was not able to run Find Surface.'
170
171     try:
172         # run the SWAF using the settings from the OVScan --&gt; check SWAF for overview experiment !!!
173         Zen.Aquisition.FindAutofocus(OVScan)
174     except:
175         print 'Was not able to run SWAF using the seetings: ', OverViewExpName
176
177 # get the resulting z-position
178 znew = Zen.Devices.Focus.ActualPosition
179
180 # try to adapt the Tile Experiment with new Z-Position
181 try:
182     Tiles.TileTools.ModifyTileRegionsZonly(OVScan, znew)
183     print 'Adapted Z-Position of Tile OverView. New Z = ',znew
184 except:
185     print 'Was not able to adapt Z-Position of Overview Scan Experiment.'
186
187 # execute the experiment
188 print 'Running OverviewScan Experiment.\n'
189 output_OVScan = Zen.Aquisition.Execute(OVScan)
190 OVScan.Close()
191 # For testing purposes - Load overview scan image automatically instead of executing the "real" experiment
192 #output_OVScan =
193 #    Zen.Application.LoadImage(r'c:\Users\MSRH\Documents\Testdata_Zeiss\RareEvent_Test_Wizard\OverViewScan_Test_raw.czii',
194 #    False)
195
196 # show the overview scan inside the document area
197 Zen.Application.Documents.Add(output_OVScan)
198 ovdoc = Zen.Application.Documents.GetByName(output_OVScan.Name)
199
200 # save the overview scan image inside the select folder
201 output_OVScan.Save(OutputFolder + '\\\\' + ovscan_name)
202
203 ##### END OVERVIEW SCAN EXPERIMENT #####
204
205 # Load analysis setting created by the wizard or an separate macro
206 ias = ZenImageAnalysisSetting()
207 # for simulation use: 000 - RareEventExample.czias
208 ias.Load(ImageAS)
209 # Analyse the image
210 Zen.Analyzing.Analyze(output_OVScan,ias)
211 # Create Zen table with results for all detected objects (parent class)
212 AllObj = Zen.Analyzing.CreateRegionsTable(output_OVScan)
213 # Create Zen table with results for each single object
214 SingleObj = Zen.Analyzing.CreateRegionTable(output_OVScan)
215
216 # check for existence of required column names
217 DetailScan = Zen.Aquisition.Experiments.GetByName(DetailExpName)
218 DetailScan.SetActive()
219 DetailIsFileExp = Tiles.TileTools.IsTilesExperiment(DetailScan)
220 out = RETools.RareEventTools.CheckColumns(SingleObj, DetailIsTileExp)
221 DetailScan.Close()
222
223 # 1st item is a bool indicating if all required columns could be found
224 columnsOK = out.Item1
225
226 if columnsOK == False:
227     print 'Execution stopped. Required Columns are missing.'
228     raise Exception('Execution stopped. Required Columns are missing.')
229
230 # 2nd item is a dictionary containing the tile properties
231 # colID: used keys = 'BCcolx' / 'BCcoly' / 'BCWidthcolx' / 'BCHeightcoly'
232 colID = out.Item2
233
234 # show and save data tables to the specified folder
235 Zen.Application.Documents.Add(AllObj)
236 Zen.Application.Documents.Add(SingleObj)
237 AllObj.Save(OutputFolder + '\\\\OverviewTable.csv')
238 SingleObj.Save(OutputFolder + '\\\\SingleObjectsTable.csv')
239
240 # check the number of detected objects = rows inside image analysis table
```



```
239 num_POI = SingleObj.RowCount
240
241 # switch to new magnification and optovar before the detailed scan starts
242 RETools.RareEventTools.SetObjectiveByName(objDT, False)
243 RETools.RareEventTools.SetOptovarByName(optoDT, False)
244
245 # determine offset between FindSurface and objects found by SWAF --&gt; required DF2
246 if CheckOffset==True:
247     # move to the first valid XY position to determine the offset
248     xpos_1st = SingleObj.GetValue(0, colID['BCcolx']) # get x stage position from list
249     ypos_1st = SingleObj.GetValue(0, colID['BCcoly']) # get y stage position from list
250     Zen.Devices.Stage.MoveTo(xpos_1st + dx_LSM, ypos_1st + dy_LSM)
251
252     # get the actual offset using FindSurface followed by SWAF --&gt; requires DF2
253     try:
254         dzFS = RETools.RareEventTools.getOffsetFindSurfaceSWAF(DetailExpName, False)
255         print 'Automatic Offset dzFS measured [micron]: ', round(dzFS, 2)
256         # store this new value inside the DF in order to use RecallFocus later on
257         print 'Store Z-Value Offset for RecallFocus.\n'
258         # store offset inside DF to use it with Recall Focus using DF2
259         Zen.Acquisition.StoreFocus()
260         useRecallFocus = True
261     except:
262         print 'Automatic Offset Check failed. Setting Offset dzFS = 0.'
263         dzFS = 0.0
264         useRecallFocus = False
265
266 elif CheckOffset==False:
267     dzFS = 0.0
268     useRecallFocus = False
269
270 if ManualOffset==True and CheckOffset==False:
271     # use the manually entered offset only when dz was not measured automatically
272     dzFS = OffsetValue
273     print 'Manual Offset [micron]: ', round(dzFS, 2)
274
275 # wait just to be sure everything is finished
276 time.sleep(2)
277
278 ##### START DETAILED SCAN EXPERIMENT #####
279
280 # execute detailed experiment at the position of every detected object
281 for i in range(0, num_POI, 1):
282
283     # get the object information from the position table
284     POI_ID = SingleObj.GetValue(i, 0) # get the ID of the object - IDs start with 2 !!!
285     xpos = SingleObj.GetValue(i, colID['BCcolx']) # get X-stage position from table
286     ypos = SingleObj.GetValue(i, colID['BCcoly']) # get Y-stage position from table
287
288     # move to the current position
289     Zen.Devices.Stage.MoveTo(xpos + dx_LSM, ypos + dy_LSM) # comment this, if one uses a simulated experiment
290     print 'Moving Stage to Object ID:', POI_ID, ' at ', round(xpos, 2), round(ypos, 2)
291
292     if useRecallFocus == False:
293         # Initial FindSurface before the Detail Scan starts
294         try:
295             Zen.Acquisition.FindSurface()
296         except:
297             'Was not able to run Find Surface.'
298             # calculate new focus position plus offset and move z-drive
299             zpos = Zen.Devices.Focus.ActualPosition + dzFS
300             print 'Move to Z-Position: ', round(zpos, 2)
301             Zen.Devices.Focus.MoveTo(zpos)
302
303     elif useRecallFocus == True:
304         # alternative solution - use RecallFocus
305         try:
306             Zen.Acquisition.RecallFocus()
307         except:
308             print 'Was not able to run Recall Focus.'
309             zpos = Zen.Devices.Focus.ActualPosition
310             print 'New z-position after Recall Focus: ', zpos
311
312     # load the predefined detailed scan experiment
313     DetailScan = Zen.Acquisition.Experiments.GetByName(DetailExpName)
314
315     # only modify the Tile Properties if required IAS features BoundWidth and BoundHeight were found
316     # if experiment is a Tile Experiment
317     if DetailIsTileExp == True:
318
319         # Modify tile center position - get bounding rectangle width & height in microns
320         bcwidth = SingleObj.GetValue(i, colID['BCWidthcolx'])
321         bcheight = SingleObj.GetValue(i, colID['BCHeightcoly'])
322         print 'Width and Height : ', str(round(bcwidth, 2)), str(round(bcheight, 2))
323         print 'Modifying Tile Properties XYZ Position and width & height.'
324         # Modify the XYZ position of the tile region on-the-fly
325         Tiles.TileTools.ModifyTileRegionsXYZ(DetailScan, xpos, ypos, zpos)
326         # Modify ConturSize for the tile according to the size of the bounding rectangle
327         Tiles.TileTools.ModifyTileRegionsSize(DetailScan, bcwidth, bcheight)
```



```
328     print 'New Tile Properties: ', round(xpos, 2), round(ypos, 2), round(zpos, 2), round(bcwidth, 2), round(bcheight, 2)
329
330     # execute the experiment
331     print 'Running Detail Scan Experiment at new XYZ position.'
332
333     # run the Detail Scan
334     output_detailscan = Zen.Acquisition.Execute(DetailScan)
335     DetailScan.Close()
336
337     # get the image data name
338     dtscan_name = output_detailscan.Name
339     # save the image data to the selected folder and close the image
340     output_detailscan.Save(OutputFolder + '\\\\' + output_detailscan.Name)
341     output_detailscan.Close()
342     # rename the CZI regarding to the object ID - Attention - IDs start with 2 !!!
343     newname_dtscan = 'DTScan_ID' + str(POI_ID) + '.czi'
344     if verbose:
345         print 'Renaming File: ' + dtscan_name + ' to: ' + newname_dtscan + '\n'
346     File.Move(OutputFolder + '\\\\' + dtscan_name, OutputFolder + '\\\\' + newname_dtscan)
347
348 ##### END DETAILED SCAN EXPERIMENT #####
349
350     # show the overview scan document again at the end
351     Zen.Application.Documents.ActiveDocument = ovdoc
352     print 'All Positions done. RareEvent Detection finished.'
```



## 6.2 Custom DLLs

The python script required two custom DLLs:

- **RETools.dll** - provides required functions for the Guided Acquisition.
- **TileTools.dll** - provides required functions to modify a ZEN Tile experiment.

Those DLLs have to be present in the ZEN program folder, which is usually located here (depending on the used software version):

"C:/Program Files/Carl Zeiss/ZEN 2/ZEN 2 (blue edition)".

## 7 ToDo's and Limitations

The current version of this tool has still limitations and room for improvement.

- No yet built-in check to avoid double detailed scans for objects close to each other within one field of view.
- Currently only modification of rectangular tile regions are supported.
- More methods to modify experiments and tile regions are planned for the next bigger ZEN Blue release.

## 8 Disclaimer

This is an application note that is free to use for everybody. Use it on your own risk.

Carl Zeiss Microscopy GmbH's ZEN software allows connection to the third party software, Python. Therefore Carl Zeiss Microscopy GmbH undertakes no warranty concerning Python, makes no representation that Python will work on your hardware, and will not be liable for any damages caused by the use of this extension. By running this example you agree to this disclaimer.